RESPIRADOR N95

Dasheng
Modelo: DTC3X

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**AJUSTE NASAL
INTERNO**

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DESCRIPCIÓN:

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LP: +591 715 66108
CBBA: +591 694 73570

SCZ: +591 691 57075
+591 773 00666
+591 691 27056

DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

NIOSH Reference: TN-14513

Mfr. Reference: SDHDTC3XAF-2

Centers for Disease Control
and Prevention (CDC)
National Institute for Occupational
Safety and Health (NIOSH)
National Personal Protective
Technology Laboratory (NPPTL)
P.O. Box 18070
Pittsburgh, PA 15236-0070
Phone: 412-386-4000
Fax: 412-386-4051
July 20, 2006

Ms. Maggie Zhong
Shanghai Dasheng Health Products Manufacture Co., Ltd.
Room 604, No. 7 Building
No.20 Handan Road
Shanghai, 200437
CHINA

Dear Ms. Zhong:

The National Institute for Occupational Safety and Health (NIOSH) has reviewed your request dated May 15, 2006. This request was for approval of the model DTC3X N95 filtering facepiece air purifying respirator. In addition, the request included the presentation of the Shanghai Dasheng Health Products Manufacture Co. Quality Manual, Edition B, dated May 10, 2003.

This request is granted. Approvals are granted only for documentation written in the English language. It is the manufacturer's responsibility to correctly translate materials desired in languages other than English. Approval number TC-84A-4329 has been assigned. The respirator is approved for protection at a N95 particulate efficiency level.

NIOSH has also reviewed the quality manual presented and finds that it meets or exceeds the minimum technical requirements for quality assurance plans outlined in Title 42, *Code of Federal Regulations* (CFR), Part 84.41(a) and on the bases of this review an approval is granted for this quality manual.

The CD enclosed with this letter contains the final respirator approval label. The abbreviated label has been accepted as submitted. The cautions and limitations, which apply to this approval, are on the approval label. Only those assemblies affected by this request, or where new approval numbers are assigned, apply to this approval action. Production approval labels cannot include information on unapproved configurations.

The approved assembly consists of the parts as listed on the approval label and the assembly matrix. Parts are to be marked with the numbers indicated on the approval label in a legible and permanent manner (marking cannot be removed without evidence of its previous presence).

This certificate of approval is not an endorsement of the respirator by NIOSH, and such endorsement shall not be stated or implied in advertisements or other publicity. However, you may publicize the fact that this respirator has met the requirements of Title 42, *Code of Federal Regulations*, Part 84.

No changes may be made to any respirators and accompanying documentation without prior written approval of NIOSH. Requests for changes must be submitted to NIOSH and a modification of this approval must be granted before changes are made.

A copy of the quality manual will be retained by NIOSH and incorporated into our files. Any future changes to this approved quality manual must be submitted to NIOSH for a modification of this approval.

Sincerely yours,

Heinz W. Ahlers
Chief, Technology Evaluation Branch
National Personal Protective Technology Laboratory

Enclosures

National Institute for Occupational Safety and Health
Technology Evaluation Branch
Certification, Evaluation, and Testing Section

TEST REPORT

Task Number: TN-14513

Prepared by: Jeremy Brannen

Date: July 14, 2006

Manufacturer: Shanghai Dasheng Health Products Manufacture Company

Item Tested: DTC3X

Background Information:

In an application dated May 15, 2006 Shanghai Dasheng Health Products Manufacture Company requested an approval of the model DTC3X N95 filtering facepiece air purifying respirator. In addition, the request included the presentation of the Shanghai Dasheng Health Products Manufacture Co. Quality Manual, Edition B, dated May 10, 2003.

Tests Assigned:

A. Exhalation Resistance Test	84.180
B. Inhalation Resistance Test	84.180
C. Sodium Chloride (NaCl) – N95	84.181

Overall Results:

The items tested passed laboratory testing.

Individual Test Results:

See accompanying test sheets.

National Institute for Occupational Safety and Health
Respirator Branch
Test Data Sheet

Task Number: TN-14513

Reference No.: CFR 84.181

Test: Sodium Chloride (NaCl) - N95

STP No.: 59

Manufacturer: Shanghai Dasheng Health Products Manufacture

Item Tested: DTC3X

Filter	Flow Rate	Initial Filter Resistance	Maximum Allowable Percent Leakage	Initial Percent Leakage	Maximum Percent Leakage	Result
1	85	8.6	5.00	0.254	0.258	PASS
2	85	9.0	5.00	0.163	0.164	PASS
3	85	9.5	5.00	0.147	0.150	PASS
4	85	9.1	5.00	0.327	0.330	PASS
5	85	9.1	5.00	0.312	0.316	PASS
6	85	9.6	5.00	0.296	0.306	PASS
7	85	9.7	5.00	0.314	0.315	PASS
8	85	9.1	5.00	0.325	0.335	PASS
9	85	9.2	5.00	0.246	0.253	PASS
10	85	9.4	5.00	0.383	0.385	PASS
11	85	9.5	5.00	0.283	0.290	PASS
12	85	10.0	5.00	0.444	0.446	PASS
13	85	9.9	5.00	0.327	0.327	PASS
14	85	9.0	5.00	0.740	0.793	PASS
15	85	8.2	5.00	0.412	0.437	PASS
16	85	9.0	5.00	0.356	0.369	PASS
17	85	8.9	5.00	0.338	0.351	PASS
18	85	8.8	5.00	0.480	0.489	PASS
19	85	9.1	5.00	0.374	0.384	PASS
20	85	8.4	5.00	0.515	0.552	PASS

Overall Result: PASS

Signature:

Engineering Technician

Date: 7/12/2006

Task Number: TN-14513

Reference No.: CFR 84.180

Test: Inhalation Resistance Test

STP No.: 7

Manufacturer: Shanghai Dasheng Health Products Manufacture

Item Tested: DTC3X

Comments:

Was all equipment verified to be in calibration throughout all testing?

☒ Yes ☐ No

Signature:

Engineering Technician

Date: 6/29/2006

National Institute for Occupational Safety and Health
Respirator Branch
Test Data Sheet

Task Number: TN-14513

Reference No.: CFR 84.180

Test: Inhalation Resistance Test

STP No.: 7

Manufacturer: Shanghai Dasheng Health Products Manufacture

Item Tested: DTC3X

Filter Type: Filter Only

Sample	Maximum Allowable Resistance (MM of H2O)	Actual Resistance (MM of H2O)	Result
	Inhalation	Inhalation	
1	35	9.4	PASS
2	35	9.1	PASS
3	35	9.1	PASS

Overall Result: PASS

Signature:

A handwritten signature in black ink, appearing to read "Jeremy B. Bannister". The signature is written in a cursive, flowing style.

Engineering Technician

Date: 6/29/2006

Task Number: TN-14513

Reference No.: CFR 84.180

Test: Exhalation Resistance Test

STP No.: 3

Manufacturer: Shanghai Dasheng Health Products Manufacture

Item Tested: DTC3X

Comments:

Was all equipment verified to be in calibration throughout all testing?

☒ Yes ☐ No

Signature:

Engineering Technician

Date: 6/29/2006

National Institute for Occupational Safety and Health
Respirator Branch
Test Data Sheet

Task Number: TN-14513

Reference No.: CFR 84.180

Test: Exhalation Resistance Test

STP No.: 3

Manufacturer: Shanghai Dasheng Health Products Manufacture

Item Tested: DTC3X

Filter Type: Filter Only

Sample	Maximum Allowable Resistance (MM of H ₂ O)	Actual Resistance (MM of H ₂ O)	Result
	Exhalation	Exhalation	
1	25	7.9	PASS
2	25	7.9	PASS
3	25	7.9	PASS

Overall Result: PASS

Signature:

A handwritten signature in black ink, appearing to read "Jeremy B. Bannister".

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SHANGHAI DASHENG HEALTH PRODUCTS MANUFACTURE CO. LTD. procedente de **CHINA**